



Website: www.ipsum.argoLims.net

At the Log In screen you may enter your Username and Password.



1-Click on COV-19 IDx

The screenshot shows the Ipsum COV-19 IDx web application interface. On the left is a sidebar with navigation links: Dashboard, Reports, Order Kit/supplies, and Contact Ipsum. The main content area has a header with the Ipsum logo and a 'Generate Req' button. Below the header is a form for entering patient and specimen information. The form fields are: COV-19 IDx (Barcode Number) (highlighted with a green circle), MRN, Location (dropdown), Date Collected (04/07/2020), Patient First Name, Patient Last Name, Date of Birth (mm/dd/yyyy), and Diagnosis Codes. To the right of the form is a large image of a biohazard bag with a barcode and a 'BIOHAZARD' label. The bag also has a 'SPECIMEN BAG' label. The bottom of the screen shows a Windows taskbar with the time 11:19 PM on 4/7/2020.

2. Enter the barcode number which is the number on our biohazard bags (scan barcode or type it.)

The screenshot shows the Ipsum COV-19 IDx web application interface. On the left is a sidebar with navigation links: Dashboard, Reports, Order Kit/supplies, and Contact Ipsum. The main content area has a header with the Ipsum logo and a 'Generate Req' button. Below the header is a form for entering patient and specimen information. The form fields are: COV-19 IDx (Barcode Number), MRN (highlighted with a green circle), Location (dropdown), Date Collected (04/07/2020), Patient First Name, Patient Last Name, Date of Birth (mm/dd/yyyy), and Diagnosis Codes. To the right of the form is a large image of a biohazard bag with a barcode and a 'BIOHAZARD' label. The bag also has a 'SPECIMEN BAG' label. Below the bag image are three sections: 'SPECIMEN COLLECTION INSTRUCTIONS' for the 2019 NOVEL CORONAVIRUS (COVID-19) TEST, 'COVID-19 TESTING USING NASOPHARYNGEAL (NP) SPECIMEN COLLECTION', and 'COVID-19 TESTING USING OROPHARYNGEAL (OP) SPECIMEN COLLECTION'. Each section has a 'Download PDF' button. The bottom of the screen shows a Windows taskbar with the time 11:19 PM on 4/7/2020.

3. Please enter your patient under investigation number.

ArgoLIMS
ipsum.argolims.net/remotecommandentry

ipsum
DIAGNOSTICS

Generate Req

Dashboard
Reports
Order Kit/supplies
Contact Ipsum

COV-19 IDx®

POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com

Copyright © 2020 Asterian, LLC - ArgoLIMS Version 4.3.2 - 1602 - API: 586

11:19 PM
4/7/2020

4. Select the location name with provider or location name and medical director.

5. The date will automatically populate for date collected.

ArgoLIMS
ipsum.argolims.net/remotecommandentry

ipsum
DIAGNOSTICS

Generate Req

Dashboard
Reports
Order Kit/supplies
Contact Ipsum

COV-19 IDx®

POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com


Copyright © 2020 Asterian, LLC - ArgoLIMS Version 4.3.2 - 1602 - API: 586

11:19 PM
4/7/2020

6. Enter First Name

7. Enter Last Name

8. Enter Date of Birth



ArgolIMS
ipsun.argolims.net/remotecommandentry

Ipsum
DIAGNOSTICS

Generate Req

COV-19 IDx (Barcode Number)
COV-19 IDx (Barcode Number)

MRN
MRN

Location
Select...

Date Collected
04/07/2020

Patient First Name
Patient First Name

Patient Last Name
Patient Last Name

Date of Birth
mm/dd/yyyy
Enter a date

Diagnosis Codes
111.22 - This is cool

Specimen Source A
Source...

SPECIMEN COLLECTION INSTRUCTIONS
2019 NOVEL CORONAVIRUS (COVID-19) TEST
Download PDF

COVID-19 TESTING USING NASOPHARYNGEAL (NP) SPECIMEN COLLECTION
Download PDF

COVID-19 TESTING USING OROPHARYNGEAL (OP) SPECIMEN COLLECTION
Download PDF


POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsundiagnostics.com

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 536

Type here to search

11:27 PM
4/7/2020

9. ICD 10 Code not required.



ArgolIMS
ipsun.argolims.net/remotecommandentry

Ipsum
DIAGNOSTICS

Generate Req

COV-19 IDx

MRN
MRN

Location
Select...

Date Collected
04/07/2020

Patient First Name
Patient First Name

Patient Last Name
Patient Last Name

Date of Birth
mm/dd/yyyy
Enter a date

Diagnosis Codes
111.22 - This is cool

Specimen Source A
Source...
Swab - Nose
Swab - Oropharyngeal

SPECIMEN COLLECTION INSTRUCTIONS
2019 NOVEL CORONAVIRUS (COVID-19) TEST
Download PDF

COVID-19 TESTING USING NASOPHARYNGEAL (NP) SPECIMEN COLLECTION
Download PDF

COVID-19 TESTING USING OROPHARYNGEAL (OP) SPECIMEN COLLECTION
Download PDF

POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsundiagnostics.com

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 536

Type here to search

11:28 PM
4/7/2020

10. Select Specimen source. There are only two options: nose or oropharyngeal.

ArgolIMS
ipsum.argolims.net/remotentry

ipsum
DIAGNOSTICS

Generate Req

Dashboard

Reports

Order Kit/supplies

Contact Ipsum

COV-19 IDx[®]

POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com

MRN
MRN

Location
Select...

Date Collected
04/07/2020

Patient First Name
Patient First Name

Patient Last Name
Patient Last Name

Date of Birth
mm/dd/yyyy
Enter a date

Diagnosis Codes
111.22 - This is cool

Specimen Source A
Source...
Required

Submit

COVID-19 TESTING USING
NASOPHARYNGEAL (NP)
SPECIMEN COLLECTION
Download PDF

COVID-19 TESTING USING
OROPHARYNGEAL (OP)
SPECIMEN COLLECTION
Download PDF

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 586

Type here to search

11:28 PM
4/7/2020

11. Click Submit

12. Click Dashboard to go to the main screen

ArgolIMS
ipsum.argolims.net

ipsum
DIAGNOSTICS

Generate Req

Dashboard

Reports

Order Kit/supplies

Contact Ipsum

COV-19 IDx[®]

POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com

All Search by Patient Name, DOB, Requisition Number, or Practice Name

Electronic Orders Paper Requisitions Filters Sort

Requisition	Last Name	First Name	DOB	MRN	Collected	Received	Reported	Report
Northstar Hospital (Smith, John DR)								
B7-18650006	Bricks	Lauren	07/18/1982	MMLBRICKS82	03/19/20			
8574956XZE	Doe	John	07/18/1975	LBRIKSO718	04/04/20			
B70010302Z	Jones	Tom	10/10/1986	MRN123	03/26/20			
B70010302X	Karen	Smith	08/28/1982	0023456	03/26/20			
XYZ123	Spanjer	Lauren	07/18/1982	MMLSMITH123	03/20/20			
B7-11112222	Tester	Joe	02/04/1980	111	03/19/20	03/19/20 11:44	03/19/20 16:11	View Report
ROBTST3	tester	test	02/04/2000	111	03/27/20			

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 586

Type here to search

11:06 PM
4/7/2020

All orders pending will show with a collected date in the column circled above.

ArgolIMS
ipsum.argolims.net

ipsum
DIAGNOSTICS

Generate Req

Search by Patient Name, DOB, Requisition Number, or Practice Name

Electronic Orders Paper Requisitions Filters Sort

Requisition	Last Name	First Name	DOB	MRN	Collected	Received	Reported	Report
Northstar Hospital (Smith, John DR)								
B7-18650006	Bricks	Lauren	07/18/1982	MMLBRICKS82	03/19/20			
B574956XZE	Doe	John	07/18/1975	LBRICKS0718	04/04/20			
B70010302Z	Jones	Tom	10/10/1986	MRN123	03/26/20			
B70010302X	Karen	Smith	08/28/1982	0023456	03/26/20			
XYZ123	Spanjer	Lauren	07/18/1982	MMLSMITH123	03/20/20			
B7-11112222	Tester	Joe	02/04/1980	111	03/19/20	03/19/20 11:44	03/19/20 16:11	View Report
ROBTST3	tester	test	02/04/2000	111	03/27/20			

POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 596

Type here to search

11:06 PM
4/7/2020

All specimens received by our lab will show a date and time as shown above.

ArgolIMS
ipsum.argolims.net

ipsum
DIAGNOSTICS

Generate Req

Search by Patient Name, DOB, Requisition Number, or Practice Name

Electronic Orders Paper Requisitions Filters Sort

Requisition	Last Name	First Name	DOB	MRN	Collected	Received	Reported	Report
Northstar Hospital (Smith, John DR)								
B7-18650006	Bricks	Lauren	07/18/1982	MMLBRICKS82	03/19/20			
B574956XZE	Doe	John	07/18/1975	LBRICKS0718	04/04/20			
B70010302Z	Jones	Tom	10/10/1986	MRN123	03/26/20			
B70010302X	Karen	Smith	08/28/1982	0023456	03/26/20			
XYZ123	Spanjer	Lauren	07/18/1982	MMLSMITH123	03/20/20			
B7-11112222	Tester	Joe	02/04/1980	111	03/19/20	03/19/20 11:44	03/19/20 16:11	View Report
ROBTST3	tester	test	02/04/2000	111	03/27/20			

POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 596

Type here to search

11:06 PM
4/7/2020

Once testing is complete, you will see a date in time here and the ability to view the report.

The screenshot shows the ArgolIMS web application interface. The search box is located at the top of the main content area, with a dropdown menu open showing search criteria: All, Patient (Name, DOB), Physician, Practice, Requisition, and Case. The table below displays patient requisitions with columns for Requisition, Last Name, First Name, DOB, MRN, Collected, Received, Reported, and Report.

Requisition	Last Name	First Name	DOB	MRN	Collected	Received	Reported	Report
Northstar Hospital (Smith, John DR)								
B7-18650006	Bricks	Lauren	07/18/1982	MMLBRICKS82	03/19/20			
B574956XZE	Doe	John	07/18/1975	LBRICKS0718	04/04/20			
B70010302Z	Jones	Tom	10/10/1986	MRN123	03/26/20			
B70010302X	Karen	Smith	08/28/1982	0023456	03/26/20			
XYZ123	Spanjer	Lauren	07/18/1982	MMLSMITH123	03/20/20			
B7-11112222	Tester	Joe	02/04/1980	111	03/19/20	03/19/20 11:44	03/19/20 16:11	View Report
ROBTST3	tester	test	02/04/2000	111	03/27/20			

The search box and criteria are located here.

The screenshot shows the ArgolIMS web application interface with the filter selection dropdown menu open. The dropdown menu is circled in green and shows options: All, Reported, Unreported, Received, and Not Received. The table below displays patient requisitions with columns for Requisition, Last Name, First Name, DOB, MRN, Collected, Received, Reported, and Report.

Requisition	Last Name	First Name	DOB	MRN	Collected	Received	Reported	Report
Northstar Hospital (Smith, John DR)								
B7-18650006	Bricks	Lauren	07/18/1982	MMLBRICKS82	03/19/20			
B574956XZE	Doe	John	07/18/1975	LBRICKS0718	04/04/20			
B70010302Z	Jones	Tom	10/10/1986	MRN123	03/26/20			
B70010302X	Karen	Smith	08/28/1982	0023456	03/26/20			
XYZ123	Spanjer	Lauren	07/18/1982	MMLSMITH123	03/20/20			
B7-11112222	Tester	Joe	02/04/1980	111	03/19/20	03/19/20 11:44	03/19/20 16:11	View Report
ROBTST3	tester	test	02/04/2000	111	03/27/20			

The filter selection gives you the ability to select reported, unreported, received and not received.

ArgolIMS
ipsum.argolims.net

ipsum
DIAGNOSTICS

Generate Req

Dashboard
Reports
Order Kit/supplies
Contact Ipsum

COV-19 IDx®

POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com

Search by Patient Name, DOB, Requisition Number, or Practice Name

Electronic Orders Paper Requisitions Filters If Sort

Requisition Last Name First Name DOB MRN Collected Received Reported Report

Northstar Hospital (Smith, John DR)

B7-18650006 Bricks Lauren 07/18/1982 MMLBRICKS82 03/19/20

8574956XZE Doe John 07/18/1975 LBRICKS0718 04/04/20

B70010302Z Jones Tom 10/10/1986 MRN123 03/26/20

B70010302X Karen Smith 08/28/1982 0023456 03/26/20

XYZ123 Spanjer Lauren 07/18/1982 MMLSMITH123 03/20/20

B7-11112222 Tester Joe 02/04/1980 111 03/19/20 03/19/20 11:44 03/19/20 16:11 View Report

ROBTST3 tester test 02/04/2000 111 03/27/20

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 586

Type here to search

11:39 PM
4/7/2020

The sort method gives you the ability to sort ascending or descending for collected, received or reported when selected.

ArgolIMS
ipsum.argolims.net

ipsum
DIAGNOSTICS

Generate Req

Dashboard
Reports
Order Kit/supplies
Contact Ipsum

COV-19 IDx®

POWERED BY
ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com

Search by Patient Name, DOB, Requisition Number, or Practice Name

Electronic Orders Paper Requisitions Filters If Sort

Requisition Last Name First Name DOB MRN Collected Received Reported Report

Northstar Hospital (Smith, John DR)

B7-18650006 Bricks Lauren 07/18/1982 MMLBRICKS82 03/19/20

8574956XZE Doe John 07/18/1975 LBRICKS0718 04/04/20

B70010302Z Jones Tom 10/10/1986 MRN123 03/26/20

B70010302X Karen Smith 08/28/1982 0023456 03/26/20

XYZ123 Spanjer Lauren 07/18/1982 MMLSMITH123 03/20/20

B7-11112222 Tester Joe 02/04/1980 111 03/19/20 03/19/20 11:44 03/19/20 16:11 View Report

ROBTST3 tester test 02/04/2000 111 03/27/20

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 586

Type here to search

11:41 PM
4/7/2020

All time search allows you to search by time frame when selected.

ArgolIMS

ipsum.argolims.net/reports

ipsum DIAGNOSTICS

Generate Req

Dashboard

Reports

Order Kit/supplies

Contact Ipsum

COV-19 IDx[®]

POWERED BY ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com

COV-19 IDx Results

Export

Start Date (Received)
2020-03-31

End Date (Received)
2020-04-07

Positive
Negative

Outcome	Collected	Received	Requisition	Location	State	Physician	Patient First Name
No rows found							

Previous Page 1 of 1 50 rows Next

positive negative

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 586

The reports tab allows you to visually see positive and negative results for all tests within a chosen time frame. You may search for your patients by requisition, location, state, physician, name, and date. You may export the report in different formats.

ArgolIMS

ipsum.argolims.net

ipsum DIAGNOSTICS

Generate Req

Dashboard

Reports

Order Kit/supplies

Contact Ipsum

COV-19 IDx[®]

POWERED BY ASTERIAN
FOR SUPPORT
877.869.2746
support@ipsumdiagnostics.com

Search by Patient Name, DOB, Requisition Number, or Practice Name

All Time

Electronic Orders Paper Requisitions

Requisition	Last Name	First Name	DOB	MRN	Collected	Received	Reported	Report
Northstar Hospital (Smith, John DR)								
B7-18650006	Bricks	Lauren	07/18/1982	MMLBRICKS82	03/19/20			
8574956KZE	Doe	John	07/18/1975	LBRICKS0718	04/04/20			
B70010302Z	Jones	Tom	10/10/1986	MRN123	03/26/20			
B70010302X	Karen	Smith	08/28/1982	0023456	03/26/20			
XYZ123	Spanjer	Lauren	07/18/1982	MMLSMITH123	03/20/20			
B7-11112222	Tester	Joe	02/04/1980	111	03/19/20	03/19/20 11:44	03/19/20 16:11	View Report
ROBTEST3	tester	test	02/04/2000	111	03/27/20			

Copyright © 2020 Asterian, LLC - ArgolIMS Version 4.3.2 - 1602 - API: 586

To order kit/supplies click the button above that will send you to our website at www.ipsumdiagnostics.com to the order form that will look like the below site or email us at supplies@ipsumdiagnostics.com.

The screenshot shows the Ipsum Diagnostics website with the URL ipsumdiagnostics.com/for-physicians/supply-order-form-covid-19/. The page features a navigation bar with links: Home, Services, Physicians, Patients, About, Contact, Login, and COV-19 IDx. The main heading is "Diagnostics" with a sub-path: Home > For Physicians > COV-19 IDx Supplies/Kits. On the left, there are buttons for "PROVIDER LOGIN" and "NEW ACCOUNT SETUP FORM", and a "More Questions?" button with a phone icon. On the right, there is a form titled "To create a new COVID-19 account, fill out the form below." with fields for "Sales Representative *" and "Practice Information" (Practice Name *).

The screenshot shows the ArgolIMS web application with the URL ipsum.argolims.net. The left sidebar contains a navigation menu with "Contact Ipsum" highlighted. The main content area displays a table of requisitions for "Northstar Hospital (Smith, John DR)". The table has columns: Requisition, Last Name, First Name, DOB, MRN, Collected, Received, Reported, and Report. A "View Report" button is visible next to the last row.

Requisition	Last Name	First Name	DOB	MRN	Collected	Received	Reported	Report
Northstar Hospital (Smith, John DR)								
B7-18650006	Bricks	Lauren	07/18/1982	MMLBRICKS82	03/19/20			
B574950XZE	Doe	John	07/18/1975	LBICKS0718	04/04/20			
B70010302Z	Jones	Tom	10/10/1986	MRN123	03/26/20			
B70010302X	Karen	Smith	08/28/1982	0023456	03/26/20			
XYZ123	Spanjer	Lauren	07/18/1982	MMLSMITH123	03/20/20			
B7-11112222	Tester	Joe	02/04/1980	111	03/19/20	03/19/20 11:44	03/19/20 16:11	View Report
ROBTST3	tester	test	02/04/2000	111	03/27/20			

If you have any questions concerning shipping for supplies contact us at 678-915-2299 or email us at shipping@ipsumdiagnostics.com.

The screenshot shows a web browser window with the URL ipsumdiagnostics.com/contact-us/. The website header includes the Ipsum Diagnostics logo and a navigation menu with links: Home, Services, Physicians, Patients, About, Contact, Login, and COV-19 IDx. The main content area is divided into two sections. On the left, a blue box titled "Contact Us" contains the following information: 8607 Roberts Drive, Suite 125, Atlanta, GA 30350; Phone: (678) 915-2299; Fax: (844) 364-8451; Email: info@ipsumdiagnostics.com; and Hours: Monday to Friday 8 am - 6 pm. On the right, a section titled "Send a Message" includes a brief instruction: "If you would like to know more about the services we offer or to get supplies for testing, please fill out the form below." Below this instruction are three input fields: "Name", "Email", and "Message". A blue "Send" button is located at the bottom right of the form.

Please be sure to let us know what frequency you would like your specimens picked up by. (Daily, Weekly or Call As Needed) You can do this by emailing Ipsum at shipping@ipsumdiagnostics.com. If you are doing a Call As Needed, please note that we utilize UPS and they will need a 1-2 hour window. Please be sure to call before 12/1PM EST. If you decide on a daily/ weekly pickup, you can decide on the time you want your specimens picked up by our courier.

When shipping your specimens, please note that there should only be 1 specimen inside each biohazard bag. Multiple biohazard bags (15-20) fit inside of the silver bag. Multiple silver bags can fit inside of the white UPS shipping bag provided to you. Please use larger ice packs if submitting larger quantities of specimens. Please let us know if you require more ice packs and we can send them to you.